MONOCRYL™ Plus Antibacterial

Poliglecaprone 25 (Monofilament) Suture

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MONOCRYL™ Plus Antibacterial

Poliglecaprone 25 (Monofilament) Suture

Sterile Synthetic Absorbable Surgical Suture

Description

MONOCRYL™ Plus Antibacterial Suture is a sterile, monofilament, synthetic, absorbable, surgical suture composed of ≥ 99.95 weight percent copolymer of glycolide and (epsilon) ϵ -caprolactone. The empirical molecular formula of the copolymer is $(C_2H_2O_2)_m (C_6H_{10}O_2)_n$. Poliglecaprone 25 copolymer has been found to be nonantigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption. MONOCRYL™ Plus Suture is available undyed and dyed with up to 0.05 weight percent D&C Violet No. 2 (Color Index 60725) to enhance visibility in the surgical field.

MONOCRYL™ Plus Suture contains Irgacare®† MP (triclosan), a broad-spectrum antibacterial agent, at no more than 2360 $\mu\text{g}/\text{m}$.

MONOCRYL™ Plus Suture is available in a range of gauge sizes and lengths, attached to needles of various types and sizes, and in presentations as described in the How Supplied section. MONOCRYL™ Plus Suture complies with the requirements of the European Pharmacopoeia (Ph. Eur.) for Sterile Synthetic Absorbable Monofilament Sutures and the United States Pharmacopoeia (USP) for Absorbable Surgical Suture, except for a slight oversize in diameter for all gauges as shown in the table below.

USP Suture Size	Metric Suture Size	Maximum Oversize
6-0	0.7	0.049 mm
5-0	1	0.033 mm
4-0	1.5	0.045 mm
3-0	2	0.067 mm
2-0	3	0.055 mm
0	3.5	0.088 mm
1	4	0.066 mm

The European Pharmacopoeia recognizes units of measure Metric and Ph. Eur. Sizes as equivalent, which is reflected on the labeling.

MONOCRYL™ Plus Suture is intended for use only by healthcare professionals who are trained in surgical suturing techniques.

Clinical benefits expected from soft tissue approximation and/or ligation are to promote wound healing by primary intention and to avoid bleeding from vessels during and after surgical procedure. A summary of safety and clinical performance can be found at the following link (upon activation):

<https://ec.europa.eu/tools/eudamed>.

Indications / Intended Use

MONOCRYL™ Plus Sutures are indicated for use in general soft tissue approximation and/or ligation, but not for use in cardiovascular or neurological tissues or ophthalmic surgery.

Contraindications

This suture (dyed and undyed), being absorbable, should not be used where extended approximation of tissues under stress is required.

Undyed MONOCRYL™ Plus Sutures, in particular, should not be used to close fascial tissue. MONOCRYL™ Plus Suture should not be used in patients with known allergic reactions (local hypersensitivity) to Irgacare® MP (triclosan).

Warnings

The safety and effectiveness of MONOCRYL™ Plus Sutures have not been established in the following areas: cardiovascular, neurological tissue and ophthalmic surgery.

Healthcare professionals should be familiar with surgical procedures and techniques involving absorbable sutures before employing MONOCRYL™ Plus Suture for wound closure, as risk of tissue separation/ wound dehiscence leading to impaired healing may vary with the site of application and the suture material used.

Healthcare professionals should consider the *in vivo* performance (under the Performance / Actions section) when selecting a suture.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

The use of MONOCRYL™ Plus Suture does not substitute normal observance of hygiene and/or otherwise needed antibiotic treatment.

As this is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the healthcare professional in the closure of sites that may undergo expansion, stretching or distension, or that may require additional support. As an absorbable suture, this suture may act transiently as a foreign body. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Like all foreign bodies, MONOCRYL™ Plus Suture may potentiate infection.

Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation that may result in device failure and/or cross-contamination that may lead to infection or transmission of blood-borne pathogens to patients and healthcare professionals.

Discard unintentionally opened / partially used / used devices and packages.

Precautions

As with any suture material, adequate knot security requires the standard surgical technique of flat, square ties with additional

throws as warranted by surgical circumstance and the experience of the healthcare professional. The use of additional throws may be particularly appropriate when knotting monofilament sutures. Improper knotting or damage to suture during use can cause extended surgery time or treatment failure and additional surgery.

Skin sutures that must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated by healthcare professional. Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with the absorption process.

Under some circumstances, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the healthcare professional.

Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur.

This suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions that may delay wound healing.

In handling this or any other suture material, care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point.

Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Healthcare professionals should exercise caution when handling surgical needles to avoid inadvertent needle stick injury that may result in transmission of blood-borne pathogens from contaminated needles. Broken needles may result in extended or additional surgeries or residual foreign bodies.

In the event of a product malfunction before use, such as a bent, broken or detached needle, or suture damage, the product should be discarded and a new one obtained to begin the procedure.

In the event of product malfunction during use, it is up to the discretion of the healthcare professional whether to continue or discontinue usage of the product and how to complete the procedure.

Adverse Reactions / Undesirable Side Effects

Adverse reactions associated with the use of this device include tissue separation/ wound dehiscence leading to impaired healing

or failure to provide adequate wound support in the closure of the sites where expansion, stretching, or distension occur, and in elderly, malnourished or debilitated patients or in patients suffering from other conditions that may delay wound healing. Other associated adverse events include calculus formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine or bile occurs, transient inflammatory tissue reaction, allergic reaction (local hypersensitivity) to Irgacare[®] MP (triclosan), and delayed absorption in tissue with poor blood supply. Transient local irritation and suture extrusion may occur at the wound site when skin sutures are left in place for greater than 7 days. Like all foreign bodies, MONOCRYL[™] Plus Suture may potentiate infection.

Suture breakage may result in hemorrhage. Suture material cutting or tearing through the tissue at the time of suture tensioning may cause soft tissue injury. Improper knotting or damage to suture during use can cause extended surgery time or treatment failure and additional surgery. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

Healthcare professionals should convey adverse reactions, undesirable side effects and risks associated with the product and the procedure to the patient and advise the patient to contact a healthcare professional in case of any deviation from the normal postoperative course. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

Magnetic Resonance Imaging (MRI) / Carcinogenic, Mutagenic, Toxic to Reproduction (CMR) / Endocrine Disrupting (ED) Safety Information

MONOCRYL[™] Plus Sutures are MR safe. No known CMR Category 1a/1b and ED substances are present at >0.1%. Category 1a/1b are defined as known or presumed human carcinogen (H340), mutagen (H350) or reproductive toxicant (H360), based on human evidence and animal studies.

Application / Instructions for Use

Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound characteristics.

Discard needles and sharps in “sharps” containers. Discard unintentionally opened / partially used / used devices and packages according to your facility’s policies and procedures concerning biohazardous materials and waste.

Performance / Actions

MONOCRYL[™] Plus Suture elicits an initial, minimal inflammatory reaction in tissues and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption occurs by means of hydrolysis, where the copolymer degrades to adipic acid that is subsequently absorbed and metabolized in the body. Absorption begins as loss of tensile strength followed by a loss of mass. Absorption is essentially complete between 91 and 119 days. Implantation studies in rats indicate that MONOCRYL[™] Plus Suture retains the percentage of original tensile strength as follows:

DYED		UNDYED	
Post Implantation	Approximate Original Strength Remaining	Post Implantation	Approximate Original Strength Remaining
7 days	60%-70%	7 days	50%-60%
14 days	30%-40%	14 days	20%-30%
All of the original tensile strength is essentially lost by 28 days post-implantation.		All of the original tensile strength is essentially lost by 21 days post-implantation.	

Using in vitro studies, MONOCRYL[™] Plus Suture has been shown to inhibit colonization of the suture by *Staphylococcus aureus*, *Staphylococcus epidermidis*, Methicillin-resistant *S. aureus*, Methicillin-resistant *S. epidermidis*, *Escherichia coli*, *Klebsiella pneumoniae* and *Enterobacter cloacae* microorganisms that are known to contribute to surgical site infections. Animal studies have demonstrated that MONOCRYL[™] Plus Suture inhibits bacterial colonization of suture after direct *in vivo* challenge with bacteria.

Sterility

MONOCRYL[™] Plus Suture is sterilized by ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged.

Storage

Store at or below 30°C (86°F). If the product is exposed outside specified storage range, contact manufacturer. Do not use after expiry date.

How Supplied

Please note that not all sizes are available in all markets. Please contact your local sales representative for size availability.

MONOCRYL[™] Plus Suture is available sterile, in sizes 6-0 through 1 (metric sizes 0.7 - 4), in a variety of lengths, with permanently attached needles.

MONOCRYL™ Plus Suture is available in one, two, or three dozen units per box.

Traceability

The following specific information can be found on the device packaging label: Catalogue number, Batch Code, expiry and manufacturing date, manufacturer name, address and the website and a Unique Device Identification bar code with the Global Trade Item Number information.

*Registered Trademark of BASF Group.

Symbols Used on Labeling

	Catalogue number
	Medical Device
	Do not use if package is damaged and consult instructions for use
	Caution
	Do not re-use
	Do not re-sterilize
	Single sterile barrier system with protective packaging inside
	Sterilized using ethylene oxide
	Batch code
	Date of manufacture
	Use-by date
	Unique Device Identifier
	Manufacturer
	Packaging unit
	Consult instructions for use or consult electronic instructions for use. EU: Call paper on demand help desk to get paper copies free of charge within 7 days.
	MR Safe

Contains a medicinal substance

Upper limit of temperature

Dyed - Absorbable - Monofilament - Uncoated - Suture

Undyed - Absorbable - Monofilament - Uncoated - Suture

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