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SARSTEDT Urin-Monovette®

Borsäure

DE	Gebrauchshinweis – SARSTEDT Urin-Monovette® Borsäure	2-7
EN	Instructions for Use – SARSTEDT Urine Monovette® Boric acid	8-13
BG	Инструкции за употреба – SARSTEDT Monovette® за урина с борна киселина	14-19
CS	Návod k obsluze – Zkumavka Urine Monovette® SARSTEDT s kyselinou boritou	20-25
DA	Brugsanvisning – SARSTEDT Urin-Monovette® borsyre	26-31
EL	Οδηγίες χρήσεως – SARSTEDT Monovette® ούρων με βορικό οξύ	32-37
ES	Instrucciones de uso – Monovette® de orina SARSTEDT con ácido bórico	38-43
ET	Kasutusjuhend – Boorhappega uriinikogumissüsteem SARSTEDT Urin-Monovette® Borsäure	44-49
FR	Mode d'emploi – Monovette® Urine avec acide borique SARSTEDT	50-55
HR	Uputa za upotrebu – SARSTEDT epruveta za urin Urin-Monovette® s bornom kiselinom	56-61
HU	Használati utasítás – SARSTEDT bórsavas vizelet Monovette®	62-67
IT	Istruzioni d'uso – SARSTEDT Monovette® urina con acido borico	68-73
KO	사용 설명서 – SARSTEDT 소변 Monovette® 봉산	74-79
LT	Naudojimo instrukcijos – SARSTEDT „Urin-Monovette®“ su boro rūgštimi	80-85
LV	Lietošanas norādes – SARSTEDT urīna Monovette® ar borskābi	86-91
NL	Gebruiksaanwijzing – SARSTEDT urine-Monovette® boorzuur	92-97
NO	Bruksanvisning – SARSTEDT Urin-Monovette® borsyre	98-103
PL	Instrukcja obsługi – Probówka do pobierania moczu Urin-Monovette® z kwasem borowym SARSTEDT	104-109
PT	Instruções de utilização – Monovette® para urina com ácido bórico SARSTEDT	110-115
RO	Instrucțiuni de utilizare – Monovette® urinară SARSTEDT cu acid boric	116-121
RU	Инструкция по применению – SARSTEDT Urin-Monovette® с борной кислотой	122-127
SK	Návod na použitie – SARSTEDT skúmavka Urin-Monovette® s kyselinou boritou	128-133
SL	Navodila za uporabo – Urinska epruveta z borovo kislino Urin-Monovette® SARSTEDT	134-139
SV	Bruksanvisning – SARSTEDT urin-Monovette® borsyra	140-145
TR	Kullanım için talimatlar – SARSTEDT İdrar-Monovette® Borik Asit	146-151

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DIN EN 14254: *In-vitro* diagnostic medical devices – Single-use receptacles for the collection of specimens, other than blood, from humans
CLSI* GP16 “Urinalysis Approved Guideline”

*CLSI (Clinical and Laboratory Standards Institute)

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European Urinalysis Guidelines: European Federation of Clinical Chemistry and Laboratory Medicine (EFCLM), „The EFCLM European Urinalysis Guideline 2023” Clinical Chemistry and Laboratory Medicine (CCLM), vol. 62, no. 9, 2024, pp. 1653-1786. <https://doi.org/10.1515/cclm-2024-0070>

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Instructions for Use – SARSTEDT Urine Monovette® Boric acid

Intended use

The Urine Monovette® Boric acid is used as a specimen container for the collection, transport, processing (e.g., by centrifugation) and storage of urine samples for microbiological *in-vitro* diagnostic tests. The Urine Monovette® Boric acid stabilises microbial growth at room temperature for up to 48 hours. The product is intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

The Urine Monovette® Boric acid consists of a transparent plastic container, a conical piston with piston rod, a green Luer screw cap with stopper and, depending on the version, an enclosed collection tip. The Urine Monovette® Boric acid is available with either a plastic or paper label and in sterile and non-sterile versions. In addition, different versions of the Urine Monovette® Boric acid are available, which are designed for different sample volumes. Some versions are compatible with the Tempus1800® transport system.

Urine Monovette® Boric acid

Name	Colour code	Description
Urine Monovette® Boric acid 3.2 ml	■	75 x 13 mm, round base
Urine Monovette® Boric acid 8.5 ml	■	92 x 15 mm, round base
Urine Monovette® Boric acid 10 ml	■	102 x 15 mm, round base

The product is prepared with boric acid. The average concentration of boric acid for a filled Urine Monovette® Boric acid is 1.5 % and stabilizes the growth of urinary pathogenic microorganisms in urine at room temperature for up to 48 hours without impairing their viability. Studies were conducted for selected microorganisms:

- *Escherichia coli*
- *Klebsiella pneumoniae*
- *Proteus mirabilis*
- *Streptococcus faecalis*
- *Candida albicans*

Information on product-specific studies is available on the SARSTEDT website:

www.sarstedt.com/literature-boricacid and www.sarstedt.com/en/products/diagnostic/urine/urine-monovette

Safety information and warnings

Note for sterile versions: Do not use the product if the packaging or the safety label is damaged.

1. General precautions: Use gloves and other general personal protective equipment to protect yourself from urine and possible exposure to pathogens transmitted by biological samples.
2. Handle all biological samples according to the guidelines and procedures in your facility. In case of direct contact with biological specimens, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases.
3. Always follow the safety guidelines and procedures of your facility.
4. The products are intended for single use. Dispose of the products and all collection utensils in containers for hazardous biological waste.
5. Do not use the products after the expiry date. The expiry date is the last day of the month and year indicated.
6. Only applies to the Urine Monovette® Boric acid: The latest version of the safety data sheets must be observed:
<https://www.sarstedt.com/safety-data-sheets>

Storage

Store the products at room temperature.

Transport

The Urine Monovette® is a primary receptacle according to the ADR (Packing instruction P650) and IATA regulations. Some versions of the Urine Monovette® are suitable for shipping with the Tempus1800® sample transportation system. These versions are listed on the respective product page on our website as suitable for the Tempus1800® transportation system.

Limitations

Significant underfilling can lead to incorrect results.

Note: Filling up to the fill level mark ensures an optimal urine-preservative ratio. A deviation of +/- 10% is permissible. It should be noted that at the start of urine collection, a few grains of the preparation may leak out, however, this does not affect the recommended minimum amount of boric acid or the stabilising effect.

The use of the stabilized urine sample for strip tests and other non-microbiological analysis methods has not been tested and may lead to false results.

Collecting and handling specimens

Read this document completely before starting the urine transfer.

Working material required for the urine transfer

1. Gloves, gown, goggles or other suitable protective wear to protect against pathogens or potentially infectious materials transmitted by urine.
2. Required number of Urine Monovettes.
3. Only for collection from an indwelling bladder catheter:
Required material for disinfecting the patient collection point on the urinary drainage system (observe facility guidelines on setting up and preparing the urine sample collection site).
4. Stand or other suitable device for placing the filled Urine Monovettes.
5. Disposal container for hazardous biological waste.

Recommended order of collection

If multiple Urine Monovettes are to be filled, the following order of collection is recommended:


Intended use	Urine Monovette®
Chemical tests, strip tests and urine particle analysis	Urine Monovette® Z (without stabiliser)
Microbiological tests	Urine Monovette® Z (without stabilizer) or Urine Monovette® Boric acid (with stabilizer)

Procedure for transferring the urine specimen to a Urine Monovette® Boric acid

A: Removal of the urine sample from a urine cup

Prepare the urine sample for withdrawal by carefully swirling the urine cup and then open the cup.


- 1 Only for individually sterile-wrapped Urine Monovettes: Open the blister package with the thumb and index finger using the peel-off technique. The transparent film must be facing upwards.
- 2 Hold the Urine Monovette® Boric acid upright with the tip pointing upwards and remove the stopper on the Luer screw cap (Fig. 2a). Set it aside for later. Place a suction tip on the Luer connection and make sure the suction tip is firmly attached (Fig. 2b).
- 3 Now dip the suction tip into the urine cup and withdraw the piston rod until the piston has reached the base line on the edge of the product label and the Urine Monovette® Boric acid is filled with urine (Fig. 3e). Take the Urine Monovette® Boric acid out of the liquid and hold it upright with the suction tip pointing upwards (Fig. 3d). To empty the suction tip, pull the piston rod downwards to the stop (Fig. 3e).
- 4 Remove the suction tip and dispose of it in a disposal container for hazardous biological waste (Fig. 4f).
- 5 Close the Urine Monovette® Boric acid securely with the stopper set aside previously. Make sure that there is no gap between the stopper and the screw cap (Fig. 5a). Only when using a Urine Monovette® Boric acid 3.2 ml: Pull the piston rod downwards until the piston audibly engages in the end position (Fig. 5b).
- 6 Break off the piston rod.
- 7 Invert the Urine Monovette® Boric acid around 10 times or until the preparation is completely dissolved. Place the Urine Monovette® in a stand.

B: Removal of the urine sample from an indwelling bladder catheter

Carry out the preparations for collecting the urine sample at the urine sample collection site following the guidelines at your facility.

- 1 Only for individually sterile wrapped Urine Monovettes:
Open the blister package with the thumb and index finger using the peel-off technique. The transparent film must be facing upwards.
- 2 Hold the Urine Monovette® Boric acid upright with the tip pointing upwards and remove the stopper on the Luer screw cap (Fig. 2a). Set it aside for later.
- 3 Penetrate the septum of the urine sample collection site and ensure the Urine Monovette® Boric acid is firmly secured with one hand (Fig. 3b).
- 4 Now, with the other hand, withdraw the piston rod back to the stop until the Urine Monovette® Boric acid is filled with urine (Fig. 4c). Take the Urine Monovette® Boric acid out of the urine sample collection site and hold it upright with the opening pointing upwards (Fig. 4d).
- 5 Close the Urine Monovette® Boric acid securely with the stopper set aside previously. Make sure that there is no gap between the stopper and the screw cap (Fig. 5a). Only when using a Urine Monovette® Boric acid 3.2 ml: Pull the piston rod downwards until the piston audibly engages in the end position (Fig. 5b).
- 6 Break off the piston rod.
- 7 Invert the Urine Monovette® Boric acid around 10 times or until the preparation is completely dissolved. Place the Urine Monovette® in a stand.

Freezing / thawing

The following general recommendations for freezing apply to Urine Monovettes:

Freezing below 0 °C

- Check beforehand whether freezing will adversely affect the urine sample or analysis (parameter stability, interference: e.g. haemolysis).

Important: Do not freeze urine samples if they are to be used for urine sediment.

Note: For the stability of the parameters, consult the instructions for use of the test reagents / analysis device manufacturers.

- The strength values of plastics are reduced in the temperature range below 0 °C. Mechanical stress should therefore be avoided.
- The conditions for freezing must be selected so that the contents of the Urine Monovette® freeze evenly or from bottom to top. The Urine Monovette® should have enough room in a holder or storage box to allow it to expand. Holders made of Styrofoam or metal are unsuitable as they can lead to expansion tears.

Freezing at -20 °C

Cool the Urine Monovette® in an upright position for a period of 45–60 min. from room temperature to +4 °C before it can be frozen at -20 °C.

Freezing below -20 °C

Deep-freezing to below -20 °C has not been tested by the manufacturer.

Due to the number of possible influencing factors, it is recommended to conduct freezing tests under routine laboratory conditions.

Thawing

Allow the Urine Monovette® to thaw for at least 45 min. at room temperature in an upright position. Mechanical stress must also be prevented in this step.

Thawing too quickly can impair the results of the analysis.

Centrifugation

Important! Centrifugation of SARSTEDT Urine Monovettes at 400 x g (RCF - relative centrifugal force - g-force) is advised.

Use up to a maximum of 2,000 x g, to be validated by the user themselves, is possible. Centrifuge inserts should be selected according to the size of the Urine Monovettes used. The relative centrifugal force has the following ratio to the selected rpm:

$$\text{RCF (g-force)} = 1.1.2 \times r \times (\text{rpm}/1000)^2$$

"RCF": "Relative centrifugal force"

"rpm": "revolutions per minute" (min⁻¹),

"r" [in cm]: "Radius of rotation from the middle of the centrifuge to the bottom of the Urine Monovette®."

Use suitable tube holders or inserts only. Centrifuging Urine Monovettes with cracks or centrifuging at an excessive centrifugal force can lead to breakage of the Urine Monovettes, which can release potentially hazardous substances.

Urine Monovettes should be centrifuged in accordance with the centrifugation conditions listed below. If other conditions are used, they must be validated by the user themselves.

It must be ensured that the Urine Monovettes fit well in the centrifuge inserts. Urine Monovettes that protrude above the insert can get caught in the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please observe the instructions for use of the centrifuge.

Caution! Do not remove broken Urine Monovettes by hand.

Information on disinfecting the centrifuge can be found in the instructions for use of the centrifuge.

Centrifugation recommendation:

The recommended RCF force (g-force) for urine specimens for sediment preparation is 400 x g for 5 minutes.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. The products must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. Dispose of potentially contaminated material according to the regulations and guidelines of the facility.
5. Please take note of the information in the safety data sheet.

Product-specific standards and guidelines as amended / updated

DIN EN 14254: *In-vitro* diagnostic medical devices - Single-use receptacles for the collection of specimens, other than blood, from humans
 CLSI* GP16 “Urinalysis Approved Guideline”

*CLSI (Clinical and Laboratory Standards Institute)

Further literature:

European Urinalysis Guidelines: European Federation of Clinical Chemistry and Laboratory Medicine (EFLM). "The EFLM European Urinalysis Guideline 2023" Clinical Chemistry and Laboratory Medicine (CCLM), vol. 62, no. 9, 2024, pp. 1653-1786. <https://doi.org/10.1515/cclm-2024-0070>

Key for symbols and labels:

	Order no.
	Batch number
	Use by
	CE marking without involvement of a notified body**
	<i>In vitro</i> diagnostic devices
	Labelling multilingualism
	Follow the instructions for use
	Risk of contamination if reused
	Keep away from sunlight
	Store in a dry place
	Do not use if packaging is damaged
	Manufacturer
	Country of manufacture

The following also applies to sterile products:

	Sterilisation by irradiation
	Do not resterilise
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

**During the transition period from the EU Directive 98/79/EC on *in vitro* diagnostic medical devices (IVDD) to the EU Regulation (EU) 2017/746 on *in vitro* diagnostic medical devices (IVDR), *in vitro* diagnostic medical devices (IVDs) of risk class A (sterile) can continue to be placed on the market as so-called legacy devices in accordance with the IVDD until 31 December 2029. These products are subject to the transitional provisions according to Article 110 of the IVDR (EU) 2017/746.

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DIN EN 14254: *In-vitro* diagnostic medical devices – Single-use receptacles for the collection of specimens, other than blood, from humans
 CLSI* GP16 “Urinalysis Approved Guideline”

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European Urinalysis Guidelines: European Federation of Clinical Chemistry and Laboratory Medicine (EFCLM), “The EFCLM European Urinalysis Guideline 2023” Clinical Chemistry and Laboratory Medicine (CCLM), vol. 62, no. 9, 2024, pp. 1653-1786. <https://doi.org/10.1515/cclm-2024-0070>

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DIN EN 14254: *In-vitro* diagnostic medical devices – Single-use receptacles for the collection of specimens, other than blood, from humans
CLSI* GP16 „Urinalysis Approved Guideline”

* CLSI (Clinical and Laboratory Standards Institute)

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European Urinalysis Guidelines: European Federation of Clinical Chemistry and Laboratory Medicine (EFLM), „The EFLM European Urinalysis Guideline 2023” Clinical Chemistry and Laboratory Medicine (CCLM), vol. 62, no. 9, 2024, pp. 1653-1786. <https://doi.org/10.1515/cclm-2024-0070>

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Norme i direktive specifične za proizvod i njihove izmjene

DIN EN 14254: *In-vitro* diagnostic medical devices - Single-use receptacles for the collection of specimens, other than blood, from humans
CLSI* GP16 "Urinalysis Approved Guideline"

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Daljnja literatura:

European Urinalysis Guidelines: European Federation of Clinical Chemistry and Laboratory Medicine (EFCLM), "The EFCLM European Urinalysis Guideline 2023" Clinical Chemistry and Laboratory Medicine (CCLM), vol. 62, no. 9, 2024, pp. 1653-1786. <https://doi.org/10.1515/cclm-2024-0070>

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DIN EN 14254: In-vitro diagnostic medical devices - Single-use receptacles for the collection of specimens, other than blood, from humans
CLSI* GP16 "Urinalysis Approved Guideline"

*CLSI(Clinical and Laboratory Standards Institute)

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European Urinalysis Guidelines: European Federation of Clinical Chemistry and Laboratory Medicine (EFLM), "The EFLM European Urinalysis Guideline 2023" Clinical Chemistry and Laboratory Medicine (CCLM), vol. 62, no. 9, 2024, pp. 1653-1786. <https://doi.org/10.1515/cclm-2024-0070>

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DIN EN 14254: *In-vitro* diagnostic medical devices - Single-use receptacles for the collection of specimens, other than blood, from humans
CLSI* GP16 “Urinalysis Approved Guideline”

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Papildoma literatūra:

European Urinalysis Guidelines: European Federation of Clinical Chemistry and Laboratory Medicine (EFCLM).. “The EFCLM European Urinalysis Guideline 2023” Clinical Chemistry and Laboratory Medicine (CCLM), vol. 62, Nr. 9, 2024, p. 1653-1786. <https://doi.org/10.1515/cclm-2024-0070>

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DIN EN 14254: *In-vitro* diagnostic medical devices - Single-use receptacles for the collection of specimens, other than blood, from humans
CLSI* "GP16 Urinalysis Approved Guideline"

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Videreførende litteratur:

European Urinalysis Guidelines: European Federation of Clinical Chemistry and Laboratory Medicine (EFCLM).. "The EFCLM European Urinalysis Guideline 2023" Clinical Chemistry and Laboratory Medicine (CCLM), vol. 62, nr. 9, 2024, s. 1653–1786. <https://doi.org/10.1515/cclm-2024-0070>

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DIN EN 14254: *In-vitro* diagnostic medical devices – Single-use receptacles for the collection of specimens, other than blood, from humans
CLSI* GP16 »Urinalysis Approved Guideline«

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Dodatna literatura:

European Urinalysis Guidelines: European Federation of Clinical Chemistry and Laboratory Medicine (EFCLM), »The EFCLM European Urinalysis Guideline 2023« Clinical Chemistry and Laboratory Medicine (CCLM), vol. 62, no. 9, 2024, pp. 1653-1786. <https://doi.org/10.1515/cclm-2024-0070>

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DIN EN 14254: *In-vitro* diagnostic medical devices – Single-use receptacles for the collection of specimens, other than blood, from humans
 CLSI* GP16 "Urinalysis Approved Guideline"

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European Urinalysis Guidelines: European Federation of Clinical Chemistry and Laboratory Medicine (EFCLM).. "The EFCLM European Urinalysis Guideline 2023" Clinical Chemistry and Laboratory Medicine (CCLM), vol. 62, no. 9, 2024, pp. 1653-1786. <https://doi.org/10.1515/cclm-2024-0070>

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