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SARSTEDT Safety-Multifly®

Kanülen / Multifly® Kanülen

DE	Gebrauchshinweis – SARSTEDT Safety-Multifly® Kanülen/Multifly® Kanülen	2 - 6
EN	Instructions for Use – SARSTEDT Safety-Multifly® Needles/Multifly® Needles	7 - 11
BG	Инструкции за употреба – SARSTEDT Safety-Multifly® Канюли/Multifly® Канюли	12 - 16
CS	Návod k obsluze – SARSTEDT Jehel Safety-Multifly®/Jehel Multifly®	17 - 21
DA	Brugsanvisning – SARSTEDT Safety-Multifly® Kanyler/Multifly® Kanyler	22 - 26
EL	Οδηγίες χρήσεως – SARSTEDT Βελόνες Safety-Multifly®/Βελόνες Multifly®	27 - 31
ES	Instrucciones de uso – SARSTEDT Aguja de seguridad Multifly®/Agujas Multifly®	32 - 36
ET	Kasutusjuhend – SARSTEDT Safety-Multifly® Kanüülid/Multifly® Kanüülid	37 - 41
FI	Käyttöohje – SARSTEDT Safety-Multifly®-Kanyylit/Multifly®-Kanyylit	42 - 46
FR	Mode d'emploi – SARSTEDT Aiguilles Safety-Multifly®/Multifly®	47 - 51
HR	Uputa za upotrebu – SARSTEDT Kanile Safety-Multifly®/Kanile Multifly®	52 - 56
HU	Használati utasítás – SARSTEDT Safety-Multifly® Kanülök/Multifly® Kanülök	57 - 61
IT	Istruzioni d'uso – SARSTEDT Aghi Safety-Multifly®/Aghi Multifly®	62 - 66
KO	사용 설명서 – SARSTEDT Safety-Multifly® 바늘/Multifly® 바늘	67 - 71
LT	Naudojimo instrukcijos – SARSTEDT „Safety-Multifly®“ Kaniulės/„Multifly®“ Kaniulės	72 - 76
LV	Lietošanas norādes – SARSTEDT Safety-Multifly® Kanīles/Multifly® Kanīles	77 - 81
NL	Gebruiksaanwijzing – SARSTEDT Safety-Multifly®-Naalden/Multifly®-Naalden	82 - 86
NO	Bruksanvisning – SARSTEDT Safety-Multifly®-Kanyler/Multifly®-Kanyler	87 - 91
PL	Instrukcja obsługi – SARSTEDT Igły Safety-Multifly®/Igły Multifly®	92 - 96
PT	Instruções de utilização – SARSTEDT Agulhas Safety-Multifly®/Agulhas Multifly®	97 - 101
RO	Instrucțiuni de utilizare – Canule SARSTEDT Safety-Multifly®/Canule SARSTEDT Multifly®	102 - 106
RU	Инструкция по применению – SARSTEDT Безопасные иглы-бабочки Safety-Multifly®/ иглы-бабочки Multifly®	107 - 111
SK	Návod na Použitie – SARSTEDT Safety-Multifly® Ihly/Multifly® Ihly	112 - 116
SL	Navodila za uporabo – SARSTEDT Kanila Safety-Multifly®/Kanila Multifly®	117 - 121
SV	Bruksanvisning – SARSTEDT Safety-Multifly®-Kanyler/Multifly®-Kanyler	122 - 126
TH	คำแนะนำในการใช้งาน – SARSTEDT หัวเข็ม Safety-Multifly®/หัวเข็ม Multifly®	127 - 131
TR	Kullanım için talimatlar – SARSTEDT Safety-Multifly® İğneler/Multifly® İğneler	132 - 136
UK	Інструкція із застосування – SARSTEDT Голки безпечні Multifly®/Голки Multifly®	137 - 141
ZH	使用说明 – SARSTEDT Safety-Multifly® 针头/Multifly® 针头	142 - 146

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Instructions for Use – SARSTEDT Safety-Multifly® Needles/ Multifly® Needles

Intended Use

The Safety-Multifly® Needle / Multifly® Needle is a sterile, single-use winged needle. The Safety-Multifly® Needle / Multifly® Needle forms a venous blood collection system when used with a SARSTEDT S-Monovette®. In addition, the Safety-Multifly® Needle / Multifly® Needle can be used for short-term gravity infusion (≤ 24 hours). The needle guard on the Safety-Multifly® Needle is used immediately after blood collection with the last S-Monovette® a short-term gravity infusion (≤ 24 hours) to cover the used needle and thus reduce the risk of accidental needle stick injuries. The product is intended for use in a professional environment and by medical professionals.

Product Description

The Safety-Multifly® Needle / Multifly® Needle is a sterile, pyrogen-free, single-use product individually packaged in a blister pack. Safety-Multifly® Needles / Multifly® Needles consist of a medical-grade $\frac{3}{4}$ inch (19 mm) long disposable stainless steel needle with a triple bevel cut, a flexible tube (80 mm or 200 mm) and an integrated needle holder (also known as a multi-adapter). A variant of the Safety-Multifly® Needles / Multifly® Needles is available without an integrated needle holder (multi-adapter) with a tube length of 80 mm or 240 mm. The products are available in needle diameters of 20, 21, 23 and 25 gauge (G). The needle guard on Safety-Multifly® Needles is used to cover the needle after use. The products contain only components that are DEHP-free and latex-free.

Combinations with other Products:

The Safety-Multifly® Needle / Multifly® Needle with integrated needle holder can be connected to the S-Monovette® by inserting the S-Monovette® into the integrated needle holder (multi-adapter) of the Safety-Multifly® Needle / Multifly® Needle and securing it with a $\frac{1}{4}$ turn clockwise. For further details, see the chapter on SAMPLING AND HANDLING or the instructions for use of the S-Monovette® (IFU 506; www.sarstedt.com/en/download/instructions-for-use).


The Safety-Multifly® Needle / Multifly® Needle with Luer connection only (= without integrated needle holder) can be connected to the S-Monovette® by first screwing the Multi-Adapter / Multi-Adapter Luer-Lock (14.1205 / 14.1205.100) into the Luer connection. Then, the S-Monovette® is inserted into the fitted Multi-Adapter / Multi-Adapter Luer-Lock and secured with a quarter turn clockwise.

The Safety-Multifly® Needle / Multifly® needle can also be connected to the Blood Gas Monovette® using the Membrane Adapter (14.1112). For further details, see the instructions for use of the Blood Gas Monovette® (IFU 510; www.sarstedt.com/en/download/instructions-for-use).

The Safety-Multifly® Needle / Multifly® Needle can be connected to an infusion system via its Luer connection. For further details, see the chapter on Sampling and Handling.


Product Variant:

Safety Multifly® Needle

Or:

Multifly® Needle

Needle with Triple Bevel

Safety Information and Warnings

1. Before use, each individual package must be checked for damage and contamination. The products must not be used if the individual packaging is damaged or contaminated. Do not use the product if the needle is bent or contaminated, or if the needle cover is missing, damaged or has fallen off.
2. The products must not be re-sterilised.
3. Follow the Instructions for Use.
4. Wear gloves and general personal protective equipment to protect yourself from transmissible pathogens.
5. Do not bend the needle.
6. Do not cover the product's safety mechanism with adhesive tape.
7. Always keep your hands behind the needle during use and disposal.
8. The needle must be kept as still as possible in the vein to avoid damaging the vein and causing pain or hematomas to the patient.
9. Do not pull forcefully on the winged needle or activate the product's safety mechanism prematurely, as this will impair its functionality.
10. Injections and infusions with a pressure infusion device are not permitted. The only permitted type of infusion is short-term gravity infusion (≤ 24 hours).
11. Before use for short-term gravity infusion (≤ 24 hours), ensure that the device is completely vented (primed).
12. Once the product's safety mechanism has been activated, it cannot be opened. Do not attempt to force it open.
13. When removing or disposing of the product, be aware of possible dripping of blood or medication. It is therefore recommended that the Safety-Multifly® Needle / Multifly® Needle be disposed of with the needle tip pointing upwards and the needle guard activated (Safety-Multifly® Needle).
14. Sharp and pointed objects (e.g. needles, adapters, lancets) must be disposed of immediately after use in approved, puncture-proof and break-proof disposal containers (e.g. Multi-Safe disposal boxes from SARSTEDT).
15. Treat all biological samples in accordance with your institution's guidelines and procedures and as potentially infectious. In the event of direct contact with biological material, seek medical advice immediately, as there is a potential risk of transmission of infectious diseases (e.g. HIV, HCV, HBV or other pathogens). Use needles with built-in needle protection to avoid needle stick injuries and thus contact with biological samples. The safety guidelines and procedures of your institution must always be followed.
16. The products are intended for single use only. Reuse is not permitted, as this poses a risk of cross-contamination of samples and infection with potentially serious pathogens for patients and users.
17. Do not spray the multi-adapter or Luer connector with disinfectant, as this may impair its functionality. The sterility of the product is guaranteed after aseptic opening of the blister pack.
18. The products must not be used after the expiry date. The expiry date is the last day of the month and year stated.

Contraindications

1. Short-term gravity infusion (≤ 24 hours) or venous blood collection must not be performed on patients with infections, burns, injuries or previous vascular surgery at the intended puncture site.
2. Medications with high venotoxicity cannot be infused via the winged needle.
3. Medicines with a strong effect on the circulatory system or medicines where the infusion rate must be controlled should not be infused via the winged needle.

General residual risks associated with venous blood collection and short-term gravity infusion (≤ 24 hours)

When applying the basic technique for blood collection and short-term gravity infusion (≤ 24 hours), the following general residual risks exist according to CLSI GP41 and WHO guidelines on drawing blood:

1. Pain
2. Bruising
3. Vasovagal Reactions
4. Infections
5. Nerve Injuries
6. Allergic reactions, primarily caused by disinfectants.

Sampling and Handling

Read these instructions for use and the instructions for use for S-Monovette® (IFU 506; www.sarstedt.com/en/download/instructions-for-use) in full before starting venepuncture.


Materials Required for Sample Collection

1. All necessary S-Monovette® blood collection tubes.
2. (Safety) Multifly® Needles (for Safety Multifly® Needles with Luer adapters, a multi-adapter (14.1205 / 14.1205.100) or a male Luer adapter for blood collection is also required)
3. Gloves, gowns, eye protection or other suitable protective clothing to protect against blood-borne pathogens or potentially infectious materials.
4. Labels for Patient Identification.
5. Disinfectant material for cleaning the collection site (follow the institution's guidelines for sterile sample collection to prepare the collection site). Do not use alcohol-based cleaning materials if the samples are to be used for blood alcohol testing.
6. Dry, Sterile Disposable Swabs.
7. Tourniquet.
8. Plasters or Bandages.
9. Sharps Container for Safe Disposal of Used Materials.

Perform venepuncture with the Safety-Multifly® Needle / Multifly® Needle.

Apply a tourniquet. Prepare the puncture site with a suitable disinfectant. Do not touch the venepuncture site after disinfection. Position the patient so that the vein is easily accessible.

1. Open the blister pack of the Safety-Multifly® Needle / Multifly® Needle aseptically using the peel-off technique with your thumb and index finger. The transparent film must face upwards. Once the disinfectant has dried, remove the protective cover from the needle, ensuring that the bevelled side of the Safety-Multifly® Needle / Multifly® Needle faces upwards.
- 2a/2b Perform venepuncture by bringing the coloured wings together at the top.
- 3 Perform venepuncture by bringing the coloured wings together upwards. Insert the S-Monovette® into the Multi-Adapter of the Safety-Multifly® Needle / Multifly® Needle (a) and secure (connect) it by turning it clockwise by a quarter (b).
- 4 Draw blood with the S-Monovette® by gently pulling the plunger rod down to the bottom of the S-Monovette®. Wait until the blood flow stops to ensure correct filling. Loosen the tourniquet during blood collection with the first S-Monovette®.
- 5 Once the blood collection is complete, first turn the S-Monovette® a quarter turn anti-clockwise (4a) and then withdraw it from the Safety-Multifly® Needle / Multifly® Needle (disconnect) (b). The Safety-Multifly® Needle / Multifly® Needle remains in the vein during blood collection and when changing the blood collection tubes. For multiple collections, subsequent S-Monovettes are connected to the Safety-Multifly® Needle / Multifly® Needles as described above, and further samples are collected. To complete the blood collection, the last S-Monovette® is disconnected from the Safety-Multifly® Needle / Multifly® Needle and step A or step B for the Safety-Multifly® Needle, as well as step C for the Multifly® Needle, is followed.


Activation of the safety mechanism of the Safety-Multify® Needle after blood collection:


A The Safety-Multify® Needle is withdrawn from the vein (**6**) and the needle guard is activated outside the vein (**7a**) by grasping the needle guard at the rear widened end with the thumb and index finger and pushing it over the needle, while the remaining fingers hold the tube as close as possible to the needle guard (**7a**). The needle must engage in the needle guard with an audible “click” (**7b**).

or

B The needle guard is activated while simultaneously withdrawing the Safety-Multify® Needle from the vein (**7a**). The needle must engage noticeably and audibly with a “click” in the needle guard (**7b**).

Activation of the needle guard – Needle outside the vein: Pull the Safety-Multify® Needle out of the vein and activate the needle guard.

Activating the needle guard – needle in the vein: Activate the needle guard while pulling out the Safety-Multify® Needle from the vein.

For Multify® Needle:

C The Multify® Needle is withdrawn from the vein (**6**).

Dispose of the Multify® Needle (**C8**, Multify® Needle pictogram) or the securely locked Safety-Multify® Needle (**A/B8**, Safety-Multify® Needle pictogram) immediately in a sharps container.

Handling for blood collection: see also the handling video for the Safety Multify® Needle: [www.sarstedt.com/en/download/videos/ Diagnostics / Venous blood](http://www.sarstedt.com/en/download/videos/Diagnostics/Venous%20blood)

Short-term gravity infusion (≤ 24 hours) with the Safety-Multify® Needle / Multify® Needle

Prepare the puncture site with a suitable disinfectant. Do not touch the venous puncture site after disinfection. Position the patient so that the vein is easily accessible. Open the blister pack of the Safety-Multify® Needle / Multify® Needle aseptically using the peel-off technique with your thumb and index finger. The transparent film must face upwards.

9 Remove the multi-adapter or cap at the end of the tube by twisting and pulling. Once the disinfectant has dried on the intended puncture site, remove the Luer plug from the infusion system and connect the female Luer connection at the end of the Safety-Multify® / Multify® Needle tube to the male Luer connection of the infusion system.

10 Pre-fill the infusion system in accordance with internal guidelines.


Maximum venting volume of the Safety-Multify® Needles / Multify® Needles with needle diameters 20G, 21G, 23G, 25G:

240 mm long tube:	0,450 ml
200 mm long tube:	0,300 ml
80 mm long tube:	0,120 ml

Remove the protective cover from the needle and ensure that the bevelled side of the Safety-Multify® Needle / Multify® Needle is facing upwards. Perform the venepuncture by bringing the coloured wings together at the top (**2a/2b**). A short-term gravity infusion (≤ 24 hours) can then be performed.

11 To end the short-term gravity infusion (≤ 24 hours), disconnect the infusion system from the Safety-Multify® Needle / Multify® Needle and then pull the Safety-Multify® Needle / Multify® Needle out of the vein by following step **A** or step **B** for the Safety-Multify® Needle and step **C** for the Multify® Needle.

Dispose of the Multify® Needle (**C8**, Multify® Needle pictogram) or the securely locked Safety-Multify® Needle (**A/B8**, Safety-Multify® Needle pictogram) immediately in a sharps container (**8**).

Storage

Store the products in a dry place at room temperature until use. Do not expose the products to direct sunlight. Exceeding the maximum recommended storage temperature may impair the quality of the product.

Disposal

1. General hygiene guidelines as well as national and local regulations for the proper disposal of contaminated or infectious materials and sharp and pointed objects must be observed and followed.
2. Personal protective equipment such as disposable gloves prevents the risk of direct contact with biological samples during disposal.
3. Contaminated products must be disposed of in sharps disposal containers.
4. Potentially contaminated consumables must be disposed of in accordance with the guidelines and regulations of the institution.

Product-specific standards and guidelines as amended/updated

CLSI* PRE02 Collection of Diagnostic Venous Blood Specimens.

*CLSI (Clinical and Laboratory Standards Institute)

WHO guidelines on drawing blood: best practice in phlebotomy

DIN 13097-4: Hypodermic needles – Part 4: Point geometry, requirements and testing.

DIN EN ISO 6009: Hypodermic needles for single use – Colour coding for identification.

DIN EN 10088-1: Stainless steels – Part 1: List of stainless steels.

DIN EN ISO 9626: Stainless steel needle tubing for the manufacture of medical devices – Requirements and test methods.

Key for Icons and Labels:

	Order no.		Manufacturer
	Batch number		Country of manufacture
	Use by		Authorised representative in Switzerland
	CE marking with the involvement of a notified body		non-pyrogenic
	Medical device		DEHP free
	Unique product identification		Applies to the United States only: "Caution: Federal law restricts this product to sale by or on the order of a licensed physician."
	Multilingual		Sterilisation with ethylene oxide gas
	Follow the instructions for use		Do not resterilise
	Risk of contamination if reused		Single sterile barrier system
	Keep away from sunlight		Single sterile barrier system with protective packaging outside
	Store in a dry place		
	Do not use if packaging is damaged		

Technical modifications reserved

All serious incidents related to the product must be reported to the manufacturer and the competent national authority.

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- CLSI* PRE02 Collection of Diagnostic Venous Blood Specimens.
- *CLSI (Clinical and Laboratory Standards Institute)
- WHO guidelines on drawing blood: best practice in phlebotomy
- DIN 13097-4: Hypodermic needles - Part 4: Point geometry, requirements and testing.
- DIN EN ISO 6009: Hypodermic needles for single use - Colour coding for identification.
- DIN EN 10088- 1: Stainless steels - Part 1: List of stainless steels.
- DIN EN ISO 9626: Stainless steel needle tubing for the manufacture of medical devices - Requirements and test methods.

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