Supa Mega Slim Cassette

EAN-0102-02A, EAN-0104-02A, EAN-0106-02A, EAN-0107-02A,
EAN-0108-02A, EAN-0109-02A, EAN-0110-02A

EMERGO EUROPE
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands

Issue Date: MAY 2023 00.984.792 V4

CellPath Ltd

80 Mochdre Enterprise Park,
Newtown, Powys, SY16 4LE, UK

T: +44 (0)1686 611333

E: sales@cellpath.com

cellpath.com

Supa Mega Slim Cassette

**EAN-0102-02A, EAN-0104-02A, EAN-0106-02A, EAN-0107-02A,
EAN-0108-02A, EAN-0109-02A, EAN-0110-02A**

A small device intended to be used in the laboratory to support / contain clinical tissue specimens (e.g., biopsy samples) to facilitate specimen processing [e.g., fixation (encapsulation in paraffin wax), dehydration, infiltration] in preparation for subsequent cytological or histological examination, and storage. It is in the form of a synthetic polymer container, with or without a lid, with drainage holes to maximize fluid exposure during tissue submersion in a parent device into which the container is loaded. This is a single-use device.

For laboratory use.

Store in a dry place. Don't store in direct sunlight.

This product does not expire under normal environmental conditions.

Marking capabilities: CellMark pens, 2B pencils.

Processing Limitations: Suitable for high temperature processing. Suitable for microwave processing.

Instructions for Use

1. Place tissue into the Supa Mega Slim Cassette. We recommend that only one sample grossed to 5mm thick is placed in the cassette.
2. Close the lid and place the cassette in a fixative ensuring it is fully submerged.
3. Process the cassette through solvents according to your usual protocol until the specimens are ready for embedding. If you are using any unusual reagents, we suggest that you carry out preliminary chemical compatibility tests before using Supa Mega Slim Cassette to contain tissue samples.
4. Carefully open the cassette.
5. Embed using the Supa Mega Slim base mould.
6. Remove the block from the mould once fully cooled and remove excess wax.

7. The Supa Mega Slim Spacer Block [JFA-0100-00A] is required to ensure the cassette is clamped in the cassette clamp adequately.

NOTE: Users of Leica Biosystems RM 2200 microtomes using the Leica Super Cassette Clamp are required to make an alteration to the base of the Leica Knife Holder if using the Supa Mega Slim base mould. If this alteration is not made only blocks cast with a Standard Supa Mega Base Mould can be sectioned.

NOTE: Users of the Eprexia rotary microtomes using the Microm Super Mega Cassette Clamp are recommended to use the Supa Mega Slim Chuck Supports [JFA-0200-63A] to ensure the cassette is adequately clamped in place.

8. Dispose of the lid in accordance with laboratory guidelines.

DO NOT REUSE THIS PRODUCT.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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