Supa Mega Cassette

EAG-0102-02A, EAG-0104-02A, EAG-0105-02A, EAG-0106-02A, EAG-0107-02A, EAG-0108-02A, EAG-0109-02A, EAG-0110-02A, EAG-0111-02A

EC REP

EMERGO EUROPE
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands

CellPath
INNOVATION IN CELLULAR PATHOLOGY

CellPath Ltd
80 Mochdre Enterprise Park,
Newtown, Powys, SY16 4LE, UK
T: +44 (0)1686 611333
E: sales@cellpath.com
cellpath.com

IVD CE

Issue Date: SEPT2023 01.537.413 V4

Supa Mega Cassette

EAG-0102-02A, EAG-0104-02A, EAG-0105-02A, EAG-0106-02A, EAG-0107-02A, EAG-0108-02A, EAG-0109-02A, EAG-0110-02A, EAG-0111-02A

A small device intended to be used in the laboratory to support / contain clinical tissue specimens (e.g., biopsy samples) to facilitate specimen processing [e.g., fixation (encapsulation in paraffin wax), dehydration, infiltration] in preparation for subsequent cytological or histological examination, and storage. It is in the form of a synthetic polymer container, with or without a lid, with drainage holes to maximize fluid exposure during tissue submersion in a parent device into which the container is loaded. This is a single-use device.

For laboratory use.

Store in a dry place. Don't store in direct sunlight.

This product does not expire under normal environmental conditions.

Marking capabilities: CellMark pens, 2B pencils.

Processing Limitations: Suitable for high temperature processing. Suitable for microwave processing.

Instructions for Use

1. Place tissue into the Supa Mega Cassette. We recommend that only one sample grossed between 5mm and 10mm thick is placed in the cassette.

2. Close the lid and place the cassette in a fixative ensuring it is fully submerged.

3. Process the cassette through solvents according to your usual protocol until the specimens are ready for embedding. If you are using any unusual reagents, we suggest that you carry out preliminary chemical compatibility tests before using Supa Mega Cassette to contain tissue samples.

4. Carefully open the cassette.

5. Embed, cool and section as per your usual protocol.

6. Dispose of the lid in accordance with laboratory guidelines.

DO NOT REUSE THIS PRODUCT.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

CellPath
INNOVATION IN CELLULAR PATHOLOGY

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