Supa Mega Mothership Cassettes

EA0-0102-02A, EA0-0104-02A, EA0-0105-02A, EA0-0106-02A, EA0-0107-02A, EA0-0108-02A, EA0-0109-02A, EA0-0110-02A

EC REP

Emergo Europe,
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands

IVD

Issue Date: JULY 2024 00.834.000 V6

CellPath
A STATLAB COMPANY

CellPath Ltd
80 Mochdre Enterprise Park,
Newtown, Powys, SY16 4LE, UK
T: +44 (0)1686 611333
E: sales@cellpath.com
cellpath.com

Supa Mega Mothership Cassettes

EAO-0102-02A, EAO-0104-02A, EAO-0105-02A, EAO-0106-02A, EAO-0107-02A, EAO-0108-02A, EAO-0109-02A, EAO-0110-02A

A small device intended to be used in the laboratory to support / contain clinical tissue specimens (e.g., biopsy samples) to facilitate specimen processing [e.g., fixation (encapsulation in paraffin wax), dehydration, infiltration] in preparation for subsequent cytological or histological examination, and storage. It is in the form of a synthetic polymer container, with or without a lid, with drainage holes to maximize fluid exposure during tissue submersion in a parent device into which the container is loaded. This is a single-use device.

For professional use

Store in a dry place. Do not store in direct sunlight. This product does not expire under normal environmental conditions.

Marking capabilities: CellMark pens, 2B pencils. All sides of the cassette can be written on.

Compatible cassettes: System II Hex, System II+, System II+ PI, System III.

Processing Limitations: Suitable for high temperature processing. Suitable for microwave processing.

NOTE: Only use with the lids provided.

Instructions for Use

1. Insert the printed ID cassette into the Supa Mega Mothership Cassette.

NOTE: The System II Hex Cassette now requires a two-step click process to complete a fully secure insertion; Click 1 – the Hex Cassette is secure but not fully inserted, Click 2 – the Hex Cassette is now fully secure and flush to the surface. All other CellPath cassettes are securely inserted on hearing one click.

2. Place the tissue into the Supa Mega Mothership Cassette. We recommend that only one sample grossed to 5mm thick is placed in the cassette.

3. Close the lid, check visually and physically to ensure that the lid is securely fitted. Place the cassette in a fixative ensuring it is fully submerged.

4. Process the cassette through solvents according to your usual protocol until the specimens are ready for embedding. If you are using any unusual reagents, we suggest that you carry out preliminary chemical compatibility tests before using System Mega Mothership Cassette to contain tissue samples.

5. Carefully open the cassette.

6. Embed using the Supa Mega Mothership Base Mould.

7. Cool and section as per your usual protocol.

NOTE: Users of Leica Biosystems RM 2200 microtomes using the Leica Super Cassette Clamp are required to make an alteration to the base of the Leica Knife Holder if using the Supa Mega Mothership Base Mould. If this alteration is not made only blocks cast with a Standard Supa Mega Base Mould can be sectioned.

8. Dispose of the lid in accordance with laboratory guidelines.

DO NOT REUSE THIS PRODUCT.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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Kompatible kassetter: System II Hex, System II+, System II+ PI, System III.

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