

# Digital Waving Rotator

360° orbital waving rotator for efficient sample mixing with accurate speed control and easy angle adjustment



The Thermo Scientific™ Digital Waving Rotator is designed for ELISA, protein synthesis, gel staining and de-staining gels, antibody staining, and other time-sensitive applications.

### Key Features

- Digital speed and time displays
- Adjustable angle: 1° to 12°
- Speed range: 5 to 120rpm
- Speed control:
  - ≤30rpm: ±1rpm
  - >30rpm: ±2rpm
- Timer range: up to 99 hr. 59 min
- Capacity: 6 x 96-well plates
- Maximum load: 9.8kg (22 lbs.)

Includes Dimpled Rubber Mat (88882102), Rubber Strips (88882103, 6/Pack).  
Use optional Rubber Mat (88882101) for plates, dishes, culture bottles, and bags.  
Warranty: 2 years | Certifications: CE




Digital Waving Rotator shown with Thermo Scientific™ Nunc™ Test Tubes on a Thermo Scientific Rubber Mat. Nunc Test Tubes sold separately.

Product	Electrical Requirements	Tray Dimensions	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Digital Waving Rotator	100 to 240V, 50/60Hz	31.6 x 27.6 x 1.9cm (12.4 x 10.9 x 0.7 in.)	38.3 x 29.5 x 22.4cm (15.1 x 11.6 x 8.8 in.)	13.5kg (29.8 lbs.)	US	<b>88882003</b>
					EU/UK/CHN	<b>88882004</b>

### Accessories

Product and Description	Quantity	Catalog Number
Rubber Mat, 32.2 x 28.2cm (12.7 x 11.1 in)	Each	<b>88882101</b>
Dimpled Rubber Mat, 32.2 x 28.2cm (12.7 x 11.1 in)	Each	<b>88882102</b>
Rubber Strips	6/pack	<b>88882103</b>
Replacement Rubber Tube	Each	<b>88882128</b>
Screw for Rubber Strips	8/Pack	<b>88882129</b>
Power Adapter, US plug (125VAC, 10A)	Each	<b>88870126</b>
Power Adapter, AUS/CHN plug (250VAC, 10A)	Each	<b>88870127</b>
Power Adapter, EU plug (250VAC, 16A)	Each	<b>88870128</b>
Power Adapter, UK plug (250VAC, 13A)	Each	<b>88870129</b>



Waving Rotator accessories shown include: Replacement Rubber Tube and Rubber Strips (upper left), Screw for Rubber Strips (upper right), Rubber Mat and Dimpled Rubber Mat (lower left and right).