



Flex 3

Depth Adjustable Safety lancet

Flex 3 Safety lancet is an all-in-one solution for emergency blood collection, featuring an adjustable design with three penetration depths in a single unit. Intended for blood glucose and other in vitro diagnostic tests, it provides healthcare professionals with a safe, reliable, and convenient tool, supporting precise and comfortable blood collection while ensuring patient safety.

Technical Specifications

Colour Code		Specification	Penetration Depth	Article No.
Blue		23G	1.30 mm - 1.80 mm - 2.30 mm	395813
Yellow		26G	1.30 mm - 1.80 mm - 2.30 mm	396813
Orange		28G	1.30 mm - 1.80 mm - 2.30 mm	397813

Shelf life 5 years

USER CONVENIENCE:
Pre-loaded and pressure activated design. The trigger position indicates the before and after use.

USER SAFETY:
Needle is fully shielded before and after use.

COMFORT:
High quality tri-bevel needle with high speed penetration.

PATIENT SAFETY:
Single-use only. Sterilised by gamma irradiation.



Transit packaging Pallet

Packaging:
Primary packaging: 100 pcs.
Transit packaging: 1,000 pcs.
Pallet: 84,000 pcs.