

Procudan



Varenummer: 140103

VARENUMMER: 140103

Dato: 15/5-2025

DATO: 16. maj 2025

Product information BULKO-TEST BOTTLE 100 ML

Range of use: Milk sample bottle. Storing test materials

Properties: Aseptic test bottle for milk
Volume 100 ml. Weight 10 g
Ø36., H 200mm

Product material: PE. Polyethylene, approved for food contact
PELD Borealis FA 5223

Producing method: Blow molding

Package size: Carton a 375 pcs, 12 carton/pallet.

Product nr: 303653

Product code: 9686 000 00

Storage: Room temperature

Country of production: Norway

Dato 15/5-2025

A handwritten signature in blue ink, appearing to read 'Torbjørn Småkasin', is positioned above the printed name.

Torbjørn Småkasin

Norsk Formblåsing AS