

Sterile latexhandsker Hounisenvarenummerserie 0860.00xx


Beskrivelse


Steril latexhandske til operationer. Pudderfri, hvilket mindsker risikoen for allergi. Den mikromatterede latex giver god fingerføling og sikkert greb. Anatomisk facon og lange ærmer med rullekant.

Tekniske data

Størrelser:	5½, 6, 6½, 7, 7½, 8, 8½, 9
Materiale:	Naturgummi-latex.
Indholdsstoffer:	ZDEC, ZDBC, ZnO, Sulfur, Anitioxidant, polymer, Titanium.
Fri for:	PVC/ftalater, thiuram og MBT.
Steril:	Elektron radiation.
Opbevaring:	Køligt (5°C til 30 °C), tørt udenfor direkte sollys.
Pakning:	Parvist.
CE-klassifikation:	Klasse II a.
Overholder normerne:	EN 455 / AQL 0,65, EN 374 (vandtæt handske med lav kemisk beskyttelse)
Kvalitetscertifikater:	EN ISO 9001:2008, EN ISO 13485:2012 + AC:2012, Medicinprodukt klasse IIa MDD 93/42/EWG. CE 0321.



 EN 374-2:2003 Specified: Level 2 (AQL < 1.5) – acc. to producer	Level	AQL
	3	< 0.65
	2	< 1.5
	1	< 4.0

 EN 374-3:2003 Test acc. to EN 374-3 40 % Sodium Hydroxide / Level 6	Level 1: > 10 min.	Level 4: > 120 min.
	Level 2: > 30 min.	Level 5: > 240 min.
	Level 3: > 60 min.	Level 6: > 480 min.

