

Multi-Safe twin plus 500

Innovative disposal box for your daily safety



A new generation for safe disposal

Certified according to DIN EN ISO 23907,
TRBA 250 and UN



Directly integrated label in 9 languages – also the option of a label as a customised solution

With innovative day closure and final closure with single-handed operation and "Press for final closure" printed on it

Ø lid: 89 mm

Height: 162 mm

Ø base: 57 mm

Filling level clearly marked due to a new label design

Easy handling and increased safety when installing the lid due to snap fittings

Integrated QR code for direct online access to the operating instructions

Also available with biohazard labelling. Item no 77.3884.105

An integrated stamp in the lid allows the insulin pen to be decoupled easily and safely upon closing


Various removal and turning aids to aid ejection



Ordering information

Order no	Description	Volume	Material	Units per case	Units per pallet	Case dimensions	Pallet dimensions
77.3884.005	Multi-Safe twin plus 500	0.5 l	PP	126	3,024	59.0 x 26.6 x 45.0 cm	120 x 80 x 195 cm
77.3884.105	Multi-Safe twin plus 500 with biohazard labelling	0.5 l	PP	126	3,024	59.0 x 26.6 x 45.0 cm	120 x 80 x 195 cm

- Individually extendible and washable brackets



Ordering information for accessories

Order no	Description	Size	Material	Units per case	Case dimensions
95.963.	Wall holder	51 x 139 mm	Stainless steel	30	38.0 x 27.8 x 14.0 cm
95.965.001	Adhesive adapter	Ø 54 mm	PP	150	25.3 x 19.3 x 19.0 cm



You can find more information on our extensive range of disposal boxes as well as our full brochure 412 at www.sarstedt.com.

SARSTEDT AG & Co. KG
P.O. Box 12 20
D-51582 Nümbrecht
Phone: +49 2293 305 0
Fax: +49 2293 305 3992
export@sarstedt.com
www.sarstedt.com