

sprout™

mini centrifuge

Ergonomic Design | Easy to Operate | Increase Efficiency

With its ergonomic design and variety of colors, the Heathrow Scientific Sprout® mini centrifuge is loaded with features that bring increased efficiency to the lab.

Tool-free snapspin™ rotors and quick-spin options make it fast and easy to use.


The lid was functionally designed to fit in the palm of the hand for easy repeat opening and quick tube extraction.

Ideal for microfilter cell separations and HPLC samples, with adapters to accommodate all major microtube and PCR tube styles (2.0, 1.5, 0.5, 0.2 mL) for customized use.


Specifications

Maximum Speed:	Approx. 6,000 rpm (fully loaded)
Maximum RCF:	2000 Xg
Power:	100 to 240 VAC 50/60 HZ 0.6A Approx. 7.3 W
Dimensions:	153 x 128 x 104 mm (5.97 x 4.99 x 4.05 in.)
One-year warranty	

4 Interchangeable Plugs: Low voltage, double insulated power adaptor



Change out rotors
with the push of a button



Accessories and Storage Case Included



HS120396



HS120395



HS120301

