




CellSafe 5™ Biopsy Cassette

Cost Effective Tissue Preparation

CellPath's patented CellSafe 5™ is a five well cassette developed for the simultaneous preparation of core biopsies, ideal for breast and prostate – two of the most costly diagnostic samples to deal with. Each well is numbered for identification and, when combined with a robust protocol, will allow fast, accurate and cost effective preparation. This reduces the pressure on laboratory staff.

Part No	Description	Colour	Quantity
EBG-0201-02A	Cellsafe 5 Biopsy Cassette	White 	100
EBG-0206-02A	Cellsafe 5 Biopsy Cassette	Blue 	100


















The cassette is available in white and blue and features an integral X-ray transparent fine mesh for more efficient processing. The blue cassette aids contrast and visualisation of light coloured biopsies, normally experienced with breast tissue.

System III Inkjet Biopsy Cassette

Specially adapted surface optimises the print quality

The System III Inkjet Biopsy Cassette comes with an attached lid with an easy to use 'break-off' hinge. Incorporating square apertures of 0.9mm², the Biopsy Cassette ensures safe accommodation of smaller biopsy specimens, whilst also enabling good levels of fluid transfer. The locking mechanism of the cassette ensures maximum tissue security during processing.

SYSTEM III INKJET BIOPSY CASSETTE (For use with Leica/Sakura Printers)				
Loose		Threaded		Colour
Part No	Quantity	Part No	Quantity	
EAI-0302-10A	500	EAI-0302-72B	2000	White 
EAI-0304-10A	500	EAI-0304-72B	2000	Green 
EAI-0305-10A	500	EAI-0305-72B	2000	Grey 
EAI-0306-10A	500	EAI-0306-72B	2000	Blue 
EAI-0307-10A	500	EAI-0307-72B	2000	Orange 
EAI-0308-10A	500	EAI-0308-72B	2000	Pink 
EAI-0309-10A	500	EAI-0309-72B	2000	Yellow 
EAI-0310-10A	500	EAI-0310-72B	2000	Red 
EAI-0311-10A	500	EAI-0311-72B	2000	Peach 
EAI-0312-10A	500	EAI-0312-72B	2000	Tan 
EAI-0315-10A	500	EAI-0315-72B	2000	Brown 
EAI-0316-10A	500	EAI-0316-72B	2000	Dark Blue 
EAI-0318-10A	500	EAI-0318-72B	2000	Dark Green 
EAI-0319-10A	500	EAI-0319-72B	2000	Lavender 
EAI-0320-10A	500	EAI-0320-72B	2000	Turquoise 



RECOMMENDED
by
CellPath
INNOVATION IN CELLULAR PATHOLOGY

Precise writing
surface optimises
print quality

The cassette has a precisely engineered writing surface which optimises the print quality on Leica and Sakura printers. It is supplied either loose or threaded for easy loading into the aforementioned cassette printers (see page 19 for packing options).

KEY BENEFITS

- 1 Compatible with both Leica and Sakura cassette printers
- 2 Optimised writing surface
- 3 Increased reagent flow
- 4 Integral fixed frangible lid
- 5 High rigidity and stability
- 6 Large colour range available