

## boso medicus exclusive

The premium blood pressure instrument - individually adjustable.



Blood pressure measurement as individual as your health: this device can be tailored more precisely than any other boso product to your personal requirements and provide numerous statistical functions for an optimal assessment of blood pressure. As easy and user friendly, as you would expect from a boso product. The measurement is performed intelligently and without additional pumping, including a motion sensor.

### DOWNLOADS:

[Leaflet blood pressure instruments](#)

[User instruction](#)

All features at one glance:

- Precise blood pressure measurement at the upper arm
- Measuring accuracy certified by the strict test reports of the European Hypertension Specialist of ESH.
- Intelligent pump-action for optimal inflation to the required pressure.
- Movement error indicator helps to avoid false readings.
- Extra-large display with brilliant contrast.
- Selectable speech function (5 languages) with volume control.
- Memory for 90 readings with date and time – Monthly overview, if 3 measurement per day | additional guest mode.
- Irregular heart beat indication.
- Display of average values for a thorough blood pressure assessment during the day (morning / afternoon).
- Evaluation of blood pressure according to WHO-classification.
- Tri-check function - automatic, three-time measurement and averaging.
- Cuff for arm circumference 22 - 42 cm.
- High-quality carrying case and batteries included.

Optionally available:

- Cuffs in various sizes:
  - XS-size | arm circumference 16-22 cm.
  - Standard size | arm circumference 22 - 32 cm.
  - XL-size | arm circumference 32 - 48 cm.
- AC-adapter for battery-free use.