

SAMU® maternity pad, non-sterile

General Product Description/Intended Purpose

Intended Purpose:

SAMU® is a single use, non-sterile hygienic maternity pad for short-term absorption of body fluids after child birth. SAMU® pad can be used by lay persons and health professionals in a clinical and home environment.

Application/Indication

SAMU® maternity pad is for absorption of body fluids after child birth.

Reference Numbers*

Packaging	P10	P20	P56
MINI CLASSIC	-	716 221	-
MIDI CLASSIC	-	-	716 212
MAXI CLASSIC	716 214	-	716 213
MAXI PLUS	716 325	-	716 313

* the full assortment may not be available in all countries

Single use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request. Reuse of the product could lead to microbiological contamination and/or mild infections of the skin.

Product disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of SAMU® should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Product performance characteristics

Characteristics	unit	MINI CLASSIC	MIDI CLASSIC	MAXI CLASSIC	MAXI PLUS
total weight:*	[g]	9	17	19.6	20.9

* ≤ ± 10% tolerance

Material characteristics

SAMU® maternity pad, non-sterile

Part	Material
fluff:	Cellulose pulp
outer nonwoven:	Hydrophilic nonwoven, polypropylene
film:	Polyethylene-film, blue (only MAXI PLUS)
thread:	Blue polypropylene ribbon (only MINI CLASSIC; MIDI CLASSIC & MAXI CLASSIC)
adhesives:	Hotmelts resp. ultrasonic welding

Product characteristics

Dimensions	unit	MINI CLASSIC	MIDI CLASSIC	MAXI CLASSIC	MAXI PLUS
total length:*	[mm]	233	300	350	350
total width:*	[mm]	90	110	110	110
length of absorbent core:*	[mm]	205	270	320	320
width of absorbent core:*	[mm]	66	108	108	108

* ≤ ± 10% tolerance

Labelling

9-digit LOT-code

e.g.: 900102822



1
year 202X

001
internal key

02
internal key

82
internal key

2
cumulated
check sum

Manufacturing date

e.g.: 

2021
year

11
month

27
day

Expiry date

e.g.: 

2026
year

11
month

27
day

Shelf life 5 years

Medical Device



Keep dry



SAMU[®] maternity pad, non-sterile

Do not reuse



CE-marking



Latest date of revision: 2021-04-06