



Tested by
Medical Departments
of Infection in Denmark



Face Shield - Hounisen

For use against virus infection

Face Shield for protection against droplet-borne sources of infection. The visor protects the face (nose, mouth and eyes) against droplets from sneezing and coughing, thereby reducing infection from other people.

The visor also helps you not to touch your face, which often happens unconsciously several times an hour and is a source of increased infection risk. The visor is an easy-to-use protective device, as it is made in one piece and is put on and discarded quickly after use. It comes in one size.

The face shield is made of virgin antifog mono layer PET, which is flexible and easy to put on. The visor can be used both with and without glasses. The forehead nobs can be turned either inside or outside for the best fit.



Item no.:	P1108.9950
Description:	Face Shield
Size:	Onesize
Material:	PET antifog mono layer
Use:	Virus protection
Approval:	CE approved
Packaging:	250 pcs. Each box

plastic and perfection
HN GROUP This face shield from Hounisen offers you a Danish-produced Face Shield made under highly controlled conditions at HN Group A/S in Billund, which is certified under ISO 13485 Medical Devices.

www.hounisen.com • Phone +45 8621 0800

