

GYNAECOLOGICAL EXAMINATION COMBI 932

ART NO: 5-932.

Our gynaecological examination chairs are designed and developed in close partnership with gynaecologists and midwives in order to come up with the very best solutions. Our aim is to make patients feel secure and provide the right ergonomics for staff at clinics and major specialist clinics alike. The height of the chair is always adjusted straight up and down, which makes working easier as the patient is always the same distance away. Custom accessories and a wide range of upholstery colours give you plenty of freedom of choice.

COMBINATION STRETCHER AND GYNAECOLOGICAL EXAMINATION CHAIR

The Gyncombi is easily converted from an examination table to a gynaecological examination chair by angling the centre section up to form a backrest, extending the patient handles and inserting the leg rests. There is no gynaecological examination equipment in the way when it is used as a table. Both the height and seat of this chair are controlled electrically, allowing you to adjust the position while still having both hands free. This chair is raised and lowered in a straight line, so the patient is always the same distance away when the chair height is adjusted.

The two backs are controlled by gas springs and the entire stretcher can be moved into shock position from both sides. The height and depth of the leg rests can be adjusted and are angled automatically to fit the patient's calves, both vertically and horizontally. The removable stainless steel rinse bowl can be disinfected in a washer disinfector. The castors make it easy to move the chair even when it is occupied by patient. Choose from four lockable or for centrally locked castors.

Leg rest covers, special leg rests and a convenient examination lamp mounted directly on the Gyncombi are available as an optional extra.



GYNAECOLOGICAL EXAMINATION COMBI 932

QUALITIES

Electric adjustable height
Electric adjustable seat
Gas spring adjustable back
Gas spring adjustable head end
Raised and lowered in a straight line
Stable and smooth base
Four lockable or four centrally locked castors
Adjustable leg rest
Patient handles, foldable
Shock position both as gynaecology chair and examination stretcher
Wide range of upholstery colours

SPECIFICATION

Patient weights	165 kg
Total width, including patient knobs	740 mm
Length, precipitated with leg rests	2115 mm
Length, precipitated without leg rest	2085 mm
Seat width	660 mm
Sitting height	520-1000 mm
Weight	92 kg

UPHOLSTERY

		
White	Beige	Dark Camel
		
Light grey	Stone	Black
		
Baby Blue	Ciel	Marine
		
Blue	Dark Green	Purple
		
Cherry	Bordeaux	

You can order free color samples

OPTIONAL EXTRA

Plastics cover for legrests/pair	5-487
Colposcope holder	5-490
Legrests Oscar Special	5-575