

EXAMINATION TABLE STANDARD, ELECTRIC

ART NO: 4-911

We offer the broadest range of examination tables available on the market; from simple standard tables to models with electric height and back adjustment and flexible adaptation based entirely on what you require. We have also developed special tables for ultrasound examinations of the heart, for example. All our examination tables move straight up and down when their height is adjusted, which means they can be placed against a wall. We use top quality castors so that our tables move easily and are easy to manoeuvre. Custom accessories and a wide range of upholstery colours give you plenty of freedom of choice.

ELECTRIC STANDARD EXAMINATION TABLE WITH EXTRA LOW HEIGHT

A quality examination table offering full stability regardless of how the patient is sitting or lying. The height is adjusted straight up and down using the electric foot pedal, which means the table can be placed against a wall. The height is well increased both in the lowest and highest position. This means easier entry and exit for the patient, while always getting a good working posture.

Examination table Standard is equipped with four lockable castors, easy to move into the shock position and with body-shaped design. Choose from 14 different colors of upholstery.

Optional extras include breathing holes and side rails.



EXAMINATION TABLE STANDARD, ELECTRIC

QUALITIES

Electric adjustable height

Extra low height

Four lockable castors

Gas spring adjustable shock position

Paper holder

Stable in all positions

Wide range of upholstery colours

SPECIFICATION

Patient weights 200 kg

Lenght 2000 mm

Widht 700 mm

Height 450-950 mm

Shock position 11°

Weight 72 kg

UPHOLSTERY



White



Beige



Dark Camel



Light grey



Stone



Black



Baby Blue



Ciel



Marine



Blue



Dark Green



Purple



Cherry



Bordeaux

You can order free color samples

OPTIONAL EXTRA

Side rails/pair, special 4-661

Breathing holes 4-657

Side rails/pair 4-663