

Anabac[®]

Autoclave deodorant



interscience

Give some air to your lab!

Laboratories are very concerned about hygiene, cleanliness and comfort of work. Try **Anabac**[®], autoclave deodorant, for a perfumed touch in your lab. **Anabac**[®] autoclave deodorant bring the discrete final touch, the refinement of cleanliness.

Classic



Well perfumed
All applications

Ref 320 100

- Mint and eucalyptus fragrance



Poma



Sweet perfume
Food industry

Ref 320 200

- Apple fragrance



Citrus



Sweet perfume
Food industry

Ref 320 300

- Lemon fragrance



Floral



Well perfumed
All applications

Ref 320 400

- Flower fragrance



Peach




Well perfumed
All applications

Ref 320 500


- Peach fragrance




Advantages



During sample sterilization, bad smell can emerge from the autoclave.



Place one capsule (or 2 if the odor is strong) in the water at the bottom of the autoclave.



After a sterilization cycle, the capsule melts, allowing the perfume to be released.

- Work comfort
- Compatible with any autoclave
- No effect on the autoclaving
- Pleasant perfume
- 100% biodegradable

Applications

- For any kind of autoclave (max: 200 L), recommended for food industry.

Specifications

- Pot of 100 capsules
- Pot weight: 160 g
- 5 fragrances: - Classic (ref: 320 100)
- Poma (ref: 320 200)
- Citrus (ref: 320 300)
- Floral (ref: 320 400)
- Peach (ref: 320 500)



Discover also our disinfection range

Applications

- Laboratory: lab equipment, glassware, benchtops, small instruments. . .
- Food production sites



Clinet® 101

Disinfecting wipes

Ideal for laboratories, hospitals, samplers, home use.

Ref 310 101



Clinet® 202

*Disinfecting wipes
in individual pack*

Ideal for laboratories, hospitals, samplers, transports, hand hygiene.

Ref 310 202



Clinet® 303

*Bactericidal and fungicidal
spray ready to use*

Ideal for large surfaces, hard to reach places, angle of walls, rough edges.

Ref 310 303



Clinet® 404 SaniFood®

*Disinfecting wipes
for food environment*

For food professionals: disinfection without rinsing of surfaces and materials in contact with food.

Ref 310 404



Clinet® Anti-algae

For laboratory water baths.

1 liter bottle with incorporated doser, Clinet® Anti-algae is used diluted at 0.1% to inhibit the growth of algae in water baths.

Ref 310 400

Distributor

interscience INTERNATIONAL

30, ch. Bois Arpents - 78860 St Nom - FRANCE
T: +33 (0)1 34 62 62 61 - F: +33 (0)1 34 62 43 03
E - mail: info@interscience.com

interscience USA & CANADA

199 Weymouth ST. - ROCKLAND - MA 02370 - USA
T: +1 781 792 2133 - F: +1 781 792 2134
E - mail: sales.usa@interscience.com

interscience CHINA

上海市徐汇区吴兴路277号锦都大厦718室-200030
电话: +86 (0)21-64739390 - 手机: +86 189 3097 0733
邮 址: sales.china@interscience.com