



Product description

Order number	83.3923.300
Product description	Cell culture plate, 48-well, PS, surface: Cell+ for sensitive adherent cells (e.g. primary cells), yellow, flat base, TC Tested (sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, DNase-/RNase-/DNA-free), 1 pc./paper/plastic packaging

Product characteristics

Label/ Print	With plastic label
Type of cap	Slip-on lid
Cell type	Challenging, adherent
Base shape	Flat base
Number of wells	48
Growth area	0,64 cm ²
Volume of work	0.5 ml
Tube	Plate
Coding colour	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Max. volume	1,27 ml
Width of product	85,2 mm
Height of product	22,4 mm
Length of product	127,8 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Cap material	Polystyrene (PS)

Purity & certification

Purity standard	TC Tested
Sterilisation	yes
Pyrogen-free	yes
Non-cytotoxic	yes
Lot or serial number	Lot number

Packaging

Minimum order qty.	50
Pcs. / inner box	1
Pcs. / box	50
Pcs. / pallet	2400
Depth of case	332 mm
Width of case	297 mm
Height of case	275 mm
Case volume	0,0271 cbm
Weight of case	3,91 kg
EAN of inner box	4038917338646
EAN of outer box	4038917338639

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

