



### **Product description**

Order number	83.3923.500
Product description	Cell culture plate, 48-well, PS, surface: hydrophobic for suspension cells, green, flat base, TC Tested (sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, DNase-/RNase-/DNA-free), 1 pc./paper/plastic packaging

### **Product characteristics**

Label/ Print	With plastic label
Type of cap	Slip-on lid
Cell type	Suspension
Base shape	Flat base
Number of wells	48
Growth area	0,64 cm <sup>2</sup>
Volume of work	0.5 ml
Tube	Plate
Coding colour	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Size**

Max. volume	1,27 ml
Width of product	85,2 mm
Height of product	22,4 mm
Length of product	127,8 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Cap material	Polystyrene (PS)

### **Purity & certification**

Purity standard	TC Tested
Sterilisation	yes
Pyrogen-free	yes
Non-cytotoxic	yes
Lot or serial number	Lot number

### **Packaging**

Minimum order qty.	50
Pcs. / inner box	1
Pcs. / box	50
Pcs. / pallet	2400
Depth of case	332 mm
Width of case	297 mm
Height of case	275 mm
Case volume	0,0271 cbm
Weight of case	3,91 kg
EAN of inner box	4038917338677
EAN of outer box	4038917338660

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

