



### **Product description**

Order number	83.3926
Product description	Cell culture plate, 96-well, PS, surface: standard TC for adherent cells, red, conical base, TC Tested (sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, DNase-/RNase-/DNA-free), 1 pc./paper/plastic packaging

### **Product characteristics**

Colour code	Red
Label/ Print	With plastic label
Cap	Lid
Cell type	Adherent
Base shape	Conical bottom
Number of wells	96
Volume of work	0,29 ml
Tube	Plate

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Size**

Width of product	85,2 mm
Length of product	127,8 mm

### **Material & colours**

Product material	Polystyrene (PS)
Cap material	Polystyrene (PS)

### **Purity & certification**

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Sterilisation	sterile
Pyrogen-free	Pyrogen-free
Lot or serial number	batched

### **Packaging**

Minimum order qty.	50
Pcs. / inner box	1
Pcs. / box	50
Pcs. / pallet	3200
Depth of case	332 mm
Width of case	297 mm
Height of case	215 mm
Case volume	0,0212 cbm
Weight of case	3,29 KG
EAN of inner box	4038917338851
EAN of outer box	4038917338844

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

