



### **Product description**

Order number	82.1581.001
Product description	Microtest plate, 96-well, flat base, polystyrene, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, DNase-/RNase-/DNA-free, labelling: shelf life and batch No., with lid, 1 pc./paper/plastic packaging

### **Product characteristics**

Label/ Print	With plastic label
Type of cap	Slip-on lid
Base shape	Flat base
Number of wells	96
Tube	Plate

### **Size**

Max. volume	390 µl
Width of product	85,2 mm
Height of product	16,55 mm
Length of product	127,8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Purity standard	sterile
Sterilisation	yes
Pyrogen-free	no
Lot or serial number	Lot number

### **Packaging**

Minimum order qty.	50
Pcs. / inner box	1
Pcs. / box	50
Pcs. / pallet	3200
Depth of case	332 mm
Width of case	297 mm
Height of case	215 mm
Case volume	0,0212 cbm
Weight of case	3,54 kg
EAN of inner box	4038917131056
EAN of outer box	4038917347952

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

